

Technical Data

Product: PIEZOVAC PV 101
Part No.: 230080V01

Vacuum connection:	G1/4 female thread / DN 16 ISO-KF adapter
Measurement range:	1200 to 0.1 mbar
Measuring units:	mbar, Torr, Pa
Display:	LCD 12 mm
Measuring cycle:	100 ms
Measuring accuracy at 1200 - 10 mbar:	0.3% of full scale
Admissible ambient temperature:	+5 to +50°C
Power supply:	9 V block battery / Plug-in mains adapter
Max. operating time:	2500 h
Protection class:	IP 40
Dimensions:	see dimension sheet
Weight:	0.2 kg

Technical data are subject to change